


MedReady User Guide

Basic Parts



HOW TO USE THE MEDREADY

1. Set the Time of Day

- a. Press and hold down the TIME SET button with one finger. Press and hold HOUR button with another finger until it advances to the correct hour.
- b. Continue to hold the TIME SET button down while you advance the minutes by pressing and holding the MIN button.
- c. Continue to hold the TIME SET button down and press the AM-PM button to select AM (no light) or PM (light).
- d. A red light in the upper right window indicates PM.



2. Set the Alarm / Dispense Times

- a. With one finger, press and hold the '1' button. With another finger, press the HOUR button, then the MIN button, then the AM/PM button until the desired time is on the screen.
- b. Look above the '1' button. If there is NOT a red light, press and hold the '1' button again for 4 seconds until the indicator light turns on.
- c. Repeat these steps for setting alarms 2-4 (if needed).
- d. Check all dispense times by pressing and releasing buttons 1-4 to make sure the times and am/pm settings are correct. If the red indicator light is not lit, the unit will not alarm or dispense.




3. Fill Tray with Medication

- Chose an empty cup to designate as the 'starter cup'. This cup will remain empty. Begin by loading the first upcoming dose into the cup to the right of the starter cup.
- Continue to the right (counterclockwise) loading the next medication doses. DO NOT skip any spaces between days.
- Do not rotate tray manually while it sits in the MedReady; you can strip the gear and damage the entire machine.



4. Place Tray into MedReady

- Gently set the tray into the medready body, centering the empty starter cup under white T-Bar.
- Align right side of empty cup with black line on rim of MedReady base.
- The white T-Bar should be in the up position.
- Place lid on unit with door in the closed position.
- Tray will rotate clockwise while dispensing.



TROUBLESHOOTING FOR THE MEDREADY

THE CLOCK DISPLAY SAYS 'OFF': The unit has been put into away mode by pressing hour button too long. Simply press and hold the 'hour' button for approximately 4 seconds and the clock should return to normal

I TOOK MY MEDICINE, BUT THE ALARM WILL NOT SHUT OFF! The white T-bar was down and/or the pill door was up when the top was put on. Unlock the blue lid and remove it. Gently slide the T-bar into the up and locked position with your finger. Replace the lid making sure the pill door is closed when you place the lid.

THE MEDREADY DOES NOT BEEP OR ROTATE AT ALL. The alarm(s) isn't turned on. Remove the blue lid and look at the 4 alarm buttons below the clock. Press and hold the one you want until the indicator light turns on.

THE ALARM RINGS, BUT THE TRAY DOES NOT ROTATE. Remove lid and press the 'Alarm Test' button. If the unit rings but does not produce any other noises or rotate, the battery is not installed correctly. Replace the lid and lock it. Flip the unit upside down and open the battery door. Connect the battery wire to the port in the corner of the battery compartment. Remove the blue lid and test again. If this doesn't help, the motor or gear may be broken and the unit needs to be replaced.

THE WRONG MEDICATION IS COMING OUT: Press the 'Alarm Test' button to ensure the tray is rotating. If it is, the medication tray was not loaded correctly. It should be loaded counterclockwise starting to the right of the empty cup in the 6:00 position of the tray.

MY MEDREADY CONNECT IS NOT CONNECTING WITH THE CALL CENTER WHEN I MISS MY MEDS. Ensure you have a MedReady connect – it should not have a flashing red light on the front. Make sure the unit has the correct time and alarm times set. Make sure the unit is not more than 100 feet away from the medical alarm it is linked too. If it still isn't working, call us at the number below for help.

An instructional video on how to set up a MedReady is available at
<http://www.mthommed.com/resources/>

If you have any questions or concerns about your MedReady, please don't hesitate to call Mountain Home Medical at 1-866-686-7504 and speak to a member of our friendly office staff on weekdays from 8:00 am – 5:30 pm.